Determination of Chloride/ Sodium Chloride in Food Samples

Titration Application M102-01

Introduction

The application describes potentiometric determination of chloride/sodium chloride ("salt") in food samples such as salt, spice mixtures, cheese, meat or tomato sauce.

Required Equipment

Apparatus

□ TL 5000/TL 7000/TL 7750

 \square 20 mL exchangeable unit (WA 20). WA 10 or WA 50 would also be suitable.

TM 235 stirrer

Electrode, Cable and Electrolyte

- AgCl 62 (item # 285102100) or AgCl 62 RG (item # 285102100) are suitable electrodes. Note: AgCl 62 features liquid electrolyte, while AgCl 62 RG features gel and cannot be refilled.
- L 1 A plug cable combination (item # 285122456).

 \Box 2 mol/L KNO₃ electrolyte (item # 285138349) *if AgCl 62 is used*.

Solutions

- **Solvent:** Distilled water
- **Titration agent:** Commercially prepared silver nitrate solution (AgNO₃) 0.1 mol/L
- **Standard:** NaCl titrimetric standard and HNO₃ 1 mol/L

Procedure

Preparation of the silver nitrate solution

It is recommended the 0.1 mol/L silver nitrate solution (AgNO $_3$) be a commercially prepared solution.

Connect electrode

Connect the electrode to the TitroLine[®] 5000, TitroLine[®] 7000, or TitroLine[®] 7750. If using the AgCl 62 electrode, *open* the refilling hole during measurement or calibration.

Standard titration

Weigh 70 to 90 mg of the NaCl standard and place in a 150 mL beaker. Add approximately 80 mL distilled water and 1 ml HNO_3 1 mol/L. Place the electrode and burette tip in the sample and start the method. The titration should stop at the equivalence point.

(W*F2) / ((EQ1-B)*M*F1)

- W = weight of the NaCl standard in grams F1 = 1
- F2 = 1000; for the conversion from milligrams to grams
- B = 0; blank value

The result is the calculated as mol/L and can be automatically transferred to the exchangeable unit WA 20.



GLP documentation

Calculation formula

Titre:	(W*F2)/((EQ1-B)*M*F1) -> WA	Mol (M):	58.44000
--------	-----------------------------	----------	----------

Sample titration

<u>Chemical Equation</u>: Ag⁺ + Cl⁻ → AgCl ↓

Calculation: (EQ1-B)*T*M*F1/(W*F2)

EQ1 = mL consumption at the equivalence point B = 0; blank value T = exact concentration of the titrant in mol/L (c * factor) M = molecular/equivalent weight of NaCl or Cl F1 = 0.1; conversion factor for % (*100/1000) W = weight of the sample in grams F2 = 1

Weigh 0.05 to 10 grams of sample according to the table below and add to a 150 mL glass beaker. Dilute the sample with distilled or deionized water up to 80 to 100 mL and add 1 mL of HNO_3 1 mol/L (or similar). Place the electrode and burette tip in the sample, adjust the stirring speed and press the "START" key. Enter sample weight and sample ID, if prompted to do so. Then push OK or ENTER to execute the titration. The titration stops at the equivalence point. The result can be read from the display, printed if a printer is connected, and exported to a USB stick.

Table 1: Sample amount using 0.1 mol/L silver nitrate solution (AgNO₃):

Salt content	Sample weight
< 0.1 %	> 10 g
0.1 - 1 %	1 – 10 g
1 - 10 %	0.1 – 2 g
10 - 50 %	0.05 - 0.1 g
50 - 100 %	0.05 g

<u>Cheese samples, butter or other solid food products:</u>

Weigh the sample in a 150 to 250 mL glass beaker and add 100 mL hot water (55 °C). For butter, boiling water is recommended. It is also recommended to use a homogenizer for better extraction of NaCl from the sample. Add 1 mL HNO₃. The warm/hot sample can be titrated directly.

<u>Result example:</u>

tration graph				
mV 250.0		Salted B	Rutter 2	
225.0				
200.0				8.36 mi 200.8 mV
175.0				1
175,0				/
150,0				
125,0				
100,0				
75,0				
50,0				
25,00	1,0 2,0	3,0 4,0	5,0 6,0 7,	0 8,0
ethod data Method name: End date:	Salt in %		Titration duratio	n: 2 m 10 s
itration data	07.03.13		End office	10/10/19
Sample ID: Start mV:	Salted Butter 2 49.7 mV		Weight: End mV:	2.93600 g 225.9 mV

GLP documentation

Salt:	(EQ1-B)*T*M*F1/(W*F2)	Mol (M):	58.44300
Blank value (B):	0.0000 ml	Titre (T):	0.10000000 (a)
Factor 1 (F1):	0.1000	Weight (W):	2.93600 g (m)
Factor 2 (F2):	1.0000	Statistics:	Off

Method Information

Method data overall view

Method name: Method type: Measured value:	Salt in % Automatic titration mV	Created at: Last modification:	03/07/13 16:30:12 03/07/13 16:40:54
Titration mode:	Dynamic	Documentation:	GLP
Dynamic:	Steep		
Measuring speed / drift:	User-defined:	minimum holding time: maximum holding time: Measuring time: Drift:	03 s 15 s 03 s 10 mV/min
Initial waiting time: Titration direction: Pretitration: End value: EQ:	0 s Increase Off Off On		-
Slope value:	User-defined	Value:	200

Dosing parameter			
Dosing speed: Maximum dosing volume:	100.00 % 50.00 ml	Filling speed:	30 s
Unit values			
Unit size: Unit ID: Reagent: Batch ID: Concentration [mol/l]: Determined at: Expire date: Opened/compounded: Test according ISO 8655: Last modification:	20ml 1296649042 AgNO3 0.1 mol/L Any Comment 0.10000 01/18/13 23:13:00 12/31/12 08/19/11 01/01/00 02/15/13 9:54:17		

Please contact our titration experts if you have any application or product questions. Thanks!

Xylem Inc. 1725 Brannum Lane Yellow Springs, Ohio 45387

Application/Technical Support: +1-845-258-1200 titration.ysi@xyleminc.com

> <u>Ordering:</u> +1-937-767-7241 orders@ysi.com



Xylem Let's Solve Water

© 2017 Xylem, Inc. All rights reserved. M102-01 0617