



Food & Beverage Series

LACTOSE MEASUREMENT IN CHEESE

Application Note 206LS
YSI Life Sciences



Life Sciences
Data for Life.™

a xylem brand

Lactose Measurement in Cheese



Life Sciences
Data for Life.™

YSI Life Sciences

Application Note 206LS

INTRODUCTION

Lactose concentrations in complex matrices such as cheese can be measured directly and quickly using the YSI Model 2900 Series Biochemistry Analyzer. Measurements are virtually unaffected by color, turbidity, density, pH, or the presence of reducing substances.

When a sample is injected into the sample chamber, the lactose diffuses into the membrane containing galactose oxidase. The lactose is immediately oxidized to hydrogen peroxide and a galactose dialdehyde derivative. The hydrogen peroxide is detected amperometrically at the platinum electrode surface. The current flow at the electrode is directly proportional to the hydrogen peroxide concentration, and hence to lactose concentration.

I. MATERIALS & SETUP

A. YSI Model 2900 Series Biochemistry Analyzer equipped with a 2702 Galactose Oxidase Membrane and 2705 Buffer.

B. Lactose standards (5.00 g/L, 25.0 g/L).

C. Connect the 2900 Series instrument to a suitable power source.

D. Perform the instrument and membrane daily checks described in the Operations Manual.

E. Volumetric glassware (Class A recommended).

F. The following instrument setup is recommended:
Sample size: 25 µL

Probe A Parameters

Chemistry	Lactose
Unit	g/L
Calibrator	5.00 g/L
End Point	30 Sec

Autocal Parameters

Temperature	1°C
Time	30 Min
Sample	5 Sam
Cal Shift	2%

II. LIMITATIONS

Because galactose and other galactosides (raffinose, stachyose) are also substrates of galactose oxidase, the presence of these in the sample will interfere with this analysis.

III. METHOD

A. Determine an appropriate dilution factor. The sample injected should contain 1 to 10 g/L lactose.

B. Accurately weigh the cheese sample to be analyzed.

Transfer the sample to a beaker and add approximately 50 mL deionized water.

Stir the sample over low heat (~55°C) for 10 to 15 minutes.

Transfer the sample to a 100 mL volumetric flask, using deionized water to aid in the complete transfer. Fill the flask to the mark with deionized water and invert to mix.

C. Calibrate the 2900 Series instrument with a 5.00 g/L lactose standard solution.

D. Check the linearity of the membrane at least once a day by injection of a lactose linearity check solution (25.0 g/L). Refer to the Operations Manual for specifications.

E. Assay the sample prepared in B by aspiration into the 2900 series instrument. The linear range of the system is 0.05 to 25.0 g/L lactose. If the value reported exceeds this, further dilution is required.

continued

F. Calibration should be done frequently as described in the Operations Manual.

IV. CALCULATIONS

To calculate % lactose, multiply the reported result by the appropriate dilution factor.

Example: 9.96 g of cheese was prepared as described. When assayed, the value reported was 2.16 g/L lactose.

$\frac{\% \text{ Lactose: } 2.16 \text{ g/L} \times 0.100}{\text{L}/9.96 \text{ g}}$	$= 0.0217 \text{ g lactose/cheese}$ $= 2.17\% \text{ (w/w)}$
--	--

ORDERING INFORMATION

YSI Part Numbers:

- 2900 Biochemistry Analyzer
- 2935 Opaque Buffer Bottle
- 2702 Galactose Oxidase Membrane Kit
- 2783 Lactose Standard Solution (5.00 g/L)
- 2784 Lactose Standard Solution (25.0 g/L)
- 2705 Buffer Kit
- 2392 NaCl Solution (for membrane installation)

For further information, please contact:

YSI Life Sciences

1725 Brannum Lane | Yellow Springs, Ohio 45387

Website: ysi.com | Email: support@ysi.com

Telephone: (937) 767-7241 | Fax: (937) 767-9320



Life Sciences
Data for Life.™

a xylem brand

YSI Life Sciences develops and manufactures scientific instruments, sensors and systems that serve a variety of scientific and industrial markets worldwide. YSI has a long history in the life sciences and bioanalytical markets, most notably with our introduction of the world's first commercial whole blood glucose analyzer in 1975. Today there are over 10,000 YSI instruments installed around the world, trusted in critical situations to provide the most accurate data in the shortest time.

ysi.com/lifesciences

YSI Life Sciences | 1725 Brannum Lane | Yellow Springs, Ohio 45387

Tel +1.800.659.8895 | +1.937.767.7241 | support@ysi.com

© 2015 Xylem, Inc. All rights reserved.

xylem
Let's Solve Water