
# 910 COD Colorimeter

INSTRUCTION MANUAL: PAGES 2 - 13

MODE D'EMPLOI: PAGES 14 - 25

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a xylem brand

# Contact Information

## **ORDERING & TECHNICAL SUPPORT**

**Telephone:** (800) 897-4151  
(937) 767-7241  
Monday through Friday, 8:00 AM to 5:00 PM ET

**Fax:** (937) 767-9353 (orders)  
(937) 767-1058 (technical support)

**Email:** [environmental@ysi.com](mailto:environmental@ysi.com)

**Mail:** YSI Incorporated  
1725 Brannum Lane  
Yellow Springs, OH 45387  
USA

**Internet:** [www.ysi.com](http://www.ysi.com)

### **When placing an order please have the following information available:**

- YSI account number (if available)
- Name and Phone Number
- Model number or brief description
- Billing and shipping address
- Quantity
- Purchase Order or Credit Card

## **SERVICE INFORMATION**

YSI has authorized service centers throughout the United States and Internationally. For the nearest service center information, please visit [www.ysi.com](http://www.ysi.com) and click 'Support' or contact YSI Technical Support directly at 800-897-4151.

When returning a product for service, include the Product Return form with cleaning certification. The form must be completely filled out for an YSI Service Center to accept the instrument for service. The Product Return form may be downloaded at [www.ysi.com](http://www.ysi.com) and clicking on the 'Support' tab.

# Table of Contents

<b>General Notes</b> .....	<b>4</b>
General Notes .....	4
Method Notes .....	4
Display Backlight .....	5
Recall of Stored Data .....	5
Replacement of batteries .....	5
<b>Methods</b> .....	<b>6</b>
COD Reagents .....	6
Preparing the Sample .....	6
Guidelines for photometric measurements .....	6
Measurement .....	7
<b>Menu options</b> .....	<b>8</b>
Menu selections .....	8
Recall of stored data .....	8
Transmitting stored data .....	9
Setting date and time .....	9
<b>Calibration Mode</b> .....	<b>10</b>
User calibration .....	10
Factory calibration reset .....	11
<b>Technical data</b> .....	<b>12</b>
Operating messages .....	13
Error codes .....	13
Battery Disposal .....	13

## **CAUTION**

The accuracy of the instrument is only valid if the instrument is used in an environment with controlled electromagnetic disturbances according to DIN 61326. Wireless devices, e.g. wireless phones, must not be used near the instrument.

# General Notes

## General Notes

- Prior to measurement ensure that the sample is suitable for analysis (no major interference) and does not require any preparation i.e. pH adjustment, filtration etc.
- Reagents are designed for use in chemical analysis only and should be kept well out of the reach of children.
- Ensure proper disposal of reagent solutions.
- Material Safety Data Sheets: available on request.

## Chemical Method Notes:

### Method:

The sample is heated at 150° C for two hours in the presence of potassium dichromate, which is a strong oxidizing agent. Organic compounds in the sample are oxidized, reducing the dichromate ion ( $\text{Cr}_2\text{O}_7^{2-}$ ) to the green chromic ion ( $\text{Cr}^{3+}$ ).

In the Low Range (LR) Test (0-150 mg/L), the amount of remaining  $\text{Cr}^{6+}$  is measured; In the Mid Range (MR) (0-1,500 mg/L) or (HR) (0-15,000 mg/L), the amount of  $\text{Cr}^{3+}$  produced is measured. The results are measured as the amount of oxygen ( $\text{O}_2$ ) consumed per liter of water sample. The LR and MR test is displayed in mg/L while the HR test shows g/L.

The COD reagent contains both mercury and silver—the mercury controls chloride interferences and the silver functions as a catalyst to speed the reaction.

### Application:

Samples can be measured if the chloride content does not exceed 1,000 mg/l (LR/MR) or 10,000 mg/l (HR).

In exceptional cases, compounds contained in the water cannot be adequately oxidized. This results in low reading, compared with the reference method.

Different methods of sampling, preparation of the sample itself and the time elapsed between taking the sample and analysis can all affect the obtained results.

## Display backlight


Press the [!] key to turn the display backlight on or off. The backlight is automatically turned off during the measurement.

## Recall of stored data

If the instrument is turned on, press the [!] key for more than 4 seconds to access the recall menu.

## Replacement of batteries:

### CAUTION:

To ensure that the instrument is waterproof:

- seal ring (E) must be in position
- battery compartment cover (B) must be fixed with the four screws

If the batteries are removed for more than one minute the date and time menu starts automatically when the photometer is turned on the next time.

# Methods

## COD Reagents

Select the appropriate vial for the desired range\*:

LR-Range:	0 – 150 mg/l, part number 251990Y, 25 count
MR-Range:	0 – 1500 mg/l, part number 251991Y, 25 count
HR-Range:	0 – 15000 mg/l, part number 251992Y, 25 count

*\*Reagents are also available in 150 count packages.*

## Preparing the Sample

Open a reaction vial with a white cap and add the specified volume.  
(Ensure that appropriate personnel safety equipment is used.)

LR-/MR:	2 ml water sample
HR:	0.2 ml water sample

Prepare a blank (Note 1) by using deionized water (TOC-free) instead of the sample  
(LR/MR: 2 ml, HR: 0.2 ml).


Replace the cap tightly. Invert the vial gently several times to mix the contents (The vial will become hot during mixing!) and digest the vials for 120 minutes in the reactor at a temperature of 150 °C.

Remove the vials from the reactor and allow them to cool down to 60 °C or less. Mix the contents by inverting each vial several times while still warm. Then allow the vials to cool to ambient temperature before measuring.

## Guidelines for photometric measurements

1. Run samples and blanks with the same batch of vials.  
The blank is stable when stored in the dark and can be used for further measurements with vials from the same batch.
2. Don't place hot vials in the adapter. Allow the vials to cool to room temperature for a minimum 45 minutes. It is recommended to leave the vials to cool over night.
3. Suspended solids in the vial lead to incorrect measurements. For this reason it is important to place the vials carefully in the adapter. The precipitant at the bottom of the sample should not be suspended.
4. Clean the outside of the vials with a towel to remove fingerprints or other marks.
5. Avoid spillage of water or reagent solution into the sample chamber because this can lead to incorrect test results.
6. Contamination of the transparent cell chamber can result in incorrect readings. Check at regular intervals and – if necessary – clean the transparent cell chamber using a moist cloth or cotton swab.
7. Large temperature differences between the instrument and the environment can lead to errors – e.g. due to the formation of condensation in the cell chamber or on the vial.
8. To avoid errors caused by stray light do not use the instrument in bright sunlight.

## Measurement

Fix the adapter for 16 mm vials on the sample chamber.

Turn the unit on using the [ON/OFF] key.

The display shows the following:

RANGE

Select the required test using the [MODE] key.

Scroll Memory (SM)

To avoid unnecessary scrolling for the required test range, the instrument memorizes the last range used before being turned off. When the instrument is turned on again, the scroll list comes up with the last used test range.

The display shows the following:

RANGE

Place the blank in the adapter (Note 1–4) making sure that the marks  $\bar{x}$  are aligned.

Blanks are specially prepared for each individual test range.

Press the [ZERO/TEST] key.

The “Range“ symbol flashes for approx. 8 seconds.

The display shows the following:

After zeroing remove the vial from the adapter.

RANGE

0.0.0

Place the sample in the adapter (Note 2–4) making sure that the marks  $\bar{x}$  are aligned.

Press the [ZERO/TEST] key.

The “Range“ symbol flashes for approx. 3 seconds.

The result appears in the display.

RANGE

RESULT

LR-/MR-range:        in mg/l

HR-range:            in g/l

The result is saved automatically.

Tolerance:  $\pm 3.5\%$  (full scale)

### **Repeating the test:**

Press the [ZERO/TEST] key again.

### **New zero calibration:**

Press the [ZERO/TEST] key for 2 seconds.

# Menu Options

## Menu selections

Press the [MODE] key and hold.

Turn the unit on using the [ON/OFF] key.  
Allow the 3 decimal points to be displayed before releasing the [MODE] key.

The [!] key allows for selection of the following menu points:

- diS      recall stored data
- Prt      printing stored data
- setting the date and time
- user calibration

The selected menu is indicated by an arrow in the display.

## diS – Recall of stored data

After confirming the selection with the [MODE] key the photometer shows the last 16 data sets in the following format (automatically proceeds every 3 seconds until result is displayed):

Number	n xx (xx: 16...1)
Year	YYYY (e.g. 2010)
Date	mm.dd (monthmonth:dayday)
Time	hh:mm (hourhour:minuteminute)
Range	Range
Result	x.xx

The [ZERO/TEST] key repeats the current data set.


The [MODE] key scrolls through all stored data sets.

Exit the menu by pressing [!] key.

PrtG

E 132

### ▲ Prt – Transmitting stored data (to Printer or PC)

Note: To print data, or to transmit to a PC, the optional **Series 900 Data Hub** (Infrared Interface Module) is required.

The **Series 900 Data Hub** and the connected printer/PC must be ready. Press the [MODE] key to start the transmitting, the instrument displays “PrtG” (Printing) for approx. 1 second followed by the number of the first data set and its transmission. All data sets will be transmitted one after the other. After finishing the instrument returns to test mode.

The print job can be cancelled by pressing the [On/Off] key. The instrument turns off.

If the instrument is not able to communicate with the IRiM, a timeout occurs after approx. 2 minutes. The error E 132 is displayed for approx. 4 seconds. Subsequently, the instrument returns to test mode (see also IRiM manual).

Mode

SET

DATE

YYYY  
(2 sec.)

Mode

Zero  
Test

!

### ▲ Setting date and time (24-hour-format)

After confirming the selection with the [MODE] key the value to be edited will be shown for 2 sec.

The setting starts with the year (YYYY) followed by the actual value to be edited. The same applies for month (mm), day (dd), hour (hh) and minutes (mm). Set the minutes first in steps of 10, press the [!] key to continue setting the minutes in steps of 1.

Increase the value by pressing the [MODE] key.

Decrease the value by pressing [ZERO/TEST] key.

Proceed to the next value to be edited by pressing [!] key.

After setting the minutes and pressing the [!] key the display will show “IS SET” and the instrument returns to the measurement mode.

# Calibration Mode

user calibration (Display in calibration mode)  
factory calibration (Display in calibration mode)

## ▼ 4 User calibration (cAL)

After confirming the selection with the [MODE] key the instrument will show CAL/LR.

Scroll through ranges using the [MODE] key.


Place the **blank** in the adapter making sure that the marks  $\times$  are aligned.

Press the [ZERO/TEST] key.

The “Range“ symbol flashes for approx. 8 seconds.

The display shows the following in alternating mode:

Perform measurement with a standard of known concentration as described.

Press the [ZERO/TEST] key.

The “Range“ symbol flashes for approx. 3 seconds.

The result is shown in the display, alternating with CAL.

If the reading corresponds with the value of the calibration standard (within the specified tolerance), exit calibration mode by pressing the [ON/OFF] key.

Changing the displayed value:


Pressing the [MODE] key once increases the displayed value by 1 digit.

Pressing the [ZERO/TEST] key once decreases the displayed value by 1 digit.

Press the corresponding key until the reading equals the value of the calibration standard.

By pressing the [ON/OFF] key, the new correction factor is calculated and stored in the user calibration software.

Confirmation of calibration (3 seconds).

**Note:** The instrument calibration in the range MR is automatically taken also for HR.

## FACTORY CALIBRATION RESET

Resetting the user calibration to the original factory calibration will reset all methods and ranges.

A user calibrated method is indicated by an arrow while the test result is displayed.

To reset the calibration press both the [MODE] and [ZERO/TEST] key and **hold**.

Turn the unit on using the [ON/OFF] key.  
Release the [MODE] and [ZERO/TEST] keys after approx. 1 second.

The following messages will appear in turn on the display:

SEL  
CAL

The factory setting is active.  
(SEL stands for Select)


or:

SEL  
cAL

Calibration has been set by the user.  
(If the user calibration is to be retained, turn the unit off using the [ON/OFF] key).

Calibration is reset to the factory setting by pressing the [MODE] key.  
The following messages will appear in turn on the display:

Mode  
SEL  
CAL  
On  
Off

Turn the unit off using the [ON/OFF] key.

# Technical Data

<b>Instrument:</b>	double wavelength, automatic wavelength selection, direct reading colorimeter
<b>Light source:</b>	LEDs, interference filters (IF) and photosensor in transparent cell chamber. Wavelength specifications of the IF: 430 nm $\Delta\lambda = 5$ nm 610 nm $\Delta\lambda = 6$ nm
<b>Wavelength accuracy:</b>	$\pm 1$ nm
<b>Photometric accuracy*:</b>	3% FS (T = 20° C – 25° C)
<b>Photometric resolution:</b>	0.01 A
<b>Power supply:</b>	4 batteries (AAA/LR 03) lifetime: approx. 17 hours or 5000 tests
<b>Auto-OFF:</b>	instrument automatically turns off 10 minutes after last keypress
<b>Display:</b>	backlit LCD (on keypress)
<b>Storage:</b>	internal ring memory for 16 data sets
<b>Serial Interface:</b>	IR interface for data transfer
<b>Time:</b>	real time clock und date
<b>Calibration:</b>	user and factory calibration resetting to factory calibration possible
<b>Dimensions:</b>	155 x 75 x 35 mm (LxWxH)
<b>Weight:</b>	approx. 260 g (incl. batteries)
<b>Ambient conditions:</b>	temperature: 5–40°C rel. humidity: 30–90% (non-condensing)
<b>Waterproof:</b>	floating; IP 68 (1 hour at 0.1 meter)
<b>CE:</b>	Certificate for Declaration of CE-Conformity

\*measured with standard solutions

To ensure maximum accuracy of test results, always use the reagent systems supplied by the instrument manufacturer.

# Operating Messages & Error Codes

Hi

Measuring range exceeded or excessive turbidity.

Lo

Result below the lowest limit of the measuring range.

Replace batteries, no further tests possible.

btLo

Battery capacity is too low for the display backlight; measurement is still possible.

Store Cal Time      Date Cal  
RESULT ▼

A user calibrated method is indicated by an arrow while the test result is displayed (see “Factory calibration reset”).

## Error codes

E27 / E28 / E29

Light absorption too great. Reasons: e.g. dirty optics.

E 10 / E 11

Calibration factor „out of range“

E 20 / E 21

Too much light reaching the detector.

E23 / E24 / E25

Too much light reaching the detector.

E 22

Battery capacity was too low during measurement. Change battery.

E 70

LR: Factory calibration incorrect / erased

E 71

LR: User calibration incorrect / erased

E 72

MR: Factory calibration incorrect / erased

E 73

MR: User calibration incorrect / erased

## Battery Disposal

The 910 COD Colorimeter is powered by alkaline batteries which the user must remove and dispose of when the batteries no longer power the instrument. Disposal requirements vary by country and region, and users are expected to understand and follow the battery disposal requirements for their specific locale.

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
# Xylem |'zīləm|

- 1) The tissue in plants that brings water upward from the roots;
- 2) a leading global water technology company.

We're 12,000 people unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.

**For more information on how Xylem can help you, go to [www.xyleminc.com](http://www.xyleminc.com)**

a xylem brand

YSI  
1725 Brannum Lane  
Yellow Springs OH 45387  
Tel +1.937.767.7241  
Fax +1.937.767.9353  
[www.ysi.com](http://www.ysi.com)  
[environmental@ysi.com](mailto:environmental@ysi.com)

910 COD Colorimeter is a trademark of Xylem Inc. or one of its subsidiaries.  
© 2012 Xylem, Inc. Version 1 10 2012