

**CRUDE FAT IN FEED – SPECIAL METHOD FOR FOAMING SAMPLES**

1. Principle

The sample is heated up with hydrochloric acid to digest protein and free bound lipids. The digestion solution is filtered and the fat remaining in the filter is extracted with petroleum ether after the drying process. The solvent is distilled and the dried residue is weighed. The fat content is calculated from the difference between the initial sample weight and the weight at the end of the analysis.

This document describes the special parameters for the hydrolysis of samples showing excessive foaming or heavy delays in boiling.

2. Method

Based on

Method book III, Feed Research, Determination of Crude Fat, Chapter 5.1.1, Governmental Method, Analysis B

3. Product Table

Feed Regulations, commented texts, Product Table 1 "Single feed" and Product Table 2 "Mixed feed"

4. Chemicals

Quality: p.a.

Water: demineralized or distilled

- 4.1. Petroleum ether 40/60
- 4.2. Hydrochloric acid HCl, c = approx. 4 mol/l
- 4.3. pH-indicator paper

5. Instruments

- 5.1. Analytica balance (0,0001g)
- 5.2. Desiccator with drying agent, e. g.: Blaugel
- 5.3. Universal Mixer e.g. Moulinex 'Moulinette'
- 5.4. HYDROTHERM
- 5.5. Extraction unit SOXTHERM micro / macro with MULTISTAT, cat. no. 13-0011 or SOXTHERM Manager, cat. no. 13-0012
- 5.6. Electric drying chamber with natural aeration and automatic control of the temperature with a precision of ± 2 °C
- 5.7. Cotton wool, chemically clean and fat-free
- 5.8. HT Weighing Paper for HYDROTHERM, cat. no. 1004939
- 5.9. Folded Filter FF240, cat. no. 1004092

6. Sample Preparation

About 3 - 5 boiling stones are put into the extraction beakers and placed into a drying chamber for about 1 hour at $100\text{ °C} \pm 3\text{ K}$. After cooling down to room temperature in the desiccator, they are weighed to 1 mg accuracy. The sample has to be conserved in order to avoid spoilage and any changes of the composition. Prior to the analysis, the sample should be at room temperature. It has to be thoroughly mixed to ensure an even distribution of the fat. The analysis should be run immediately after the mixing.

7. Procedure**7.1. Hydrolysis**

The weight of the sample should be determined by aiming to achieve a target of extracting 0.5 – 1.5 g fat; so up to 10 g sample is weighed into a weighing paper (5.8.) and put into the digestion beaker (± 0.1 mg precision).

The digestion beaker containing the sample is inserted in the HYDROTHERM and is locked.

For further information regarding the sample weight please see section 9 (Product table).

Start use of the HYDROTHERM unit following the instruction manual.

A dry folded filter (5.9.) is placed into the respective position of the instrument. Then the apparatus is closed and the program can be started.



CRUDE FAT IN FEED – SPECIAL METHOD FOR FOAMING SAMPLES

The hydrochloric acid (4.2.) is added automatically. The liquid is quickly brought to boil and simmered with reduced heating capacity for about 1 h. At the end of the hydrolysis the digestion mixture is diluted with hot water to the double amount and is then immediately filtered through a pleated filter, which has been moistened automatically by the system with water. The beaker, the condenser and the filtration components are rinsed several times with hot water. The filter is rinsed with hot water till the backwash water has a neutral reaction. HYDROTHERM performs all these steps automatically.

Recommended Settings and Parameters: HYDROTHERM Soft Method

| Parameter | Setting | Unit | Comment |
|------------------------------------|---------|------|--------------------------|
| Fill levels | | | |
| HCl amount | 100 | ml | |
| H ₂ O amount / dilution | 100 | ml | |
| Heating / cooling phases | | | |
| Heat up phase | 6 | min | |
| Boiling phase I | 20 | min | Boiling phase power 50 % |
| Boiling phase II | 30 | min | Boiling phase power 30 % |
| Cool down phase - duration | 15 | min | |
| Filter moisture | | | |
| Number of moisture cycles | 3 | | |
| Moisture amount per cycle | 30 | ml | |
| Filter phase: | | | |
| Filter wait time | 5 | sec | |
| Rinsing cycles | 16 | | |
| Pipe opening time | 200 | ms | |
| Sample rinse time | 10 | sec | |
| Sample shower - amount | 40 | ml | |
| Cooling shower - amount | 20 | ml | |
| Filter shower - amount | 40 | ml | |

After the program has finished, the filter is placed on a watch glass and dried for up to 1.5 h at 103 °C ± 2 K in a drying oven. Since humidity may have an impact on the result, the drying time might have to be prolonged.

7.2. Extraction

After the cooling period, the filter is put into an extraction thimble and covered with cotton wool (5.7.). Any remaining fat traces on the watch glass have to be taken up with the cotton (5.7.), soaked with the extraction agent, and put into the extraction thimble as well. After adding 150 ml of extraction agent in case of SOXTHERM macro and 100 ml of extraction agent (4.1.) respectively in case of SOXTHERM micro, the sample is extracted using the following program:

- Solvent: Petroleum ether
- Boiling Range: 40 to 60 °C
- Amount of Solvent: 100 ml SOXTHERM micro / 150 ml SOXTHERM macro

- Size of Extraction Beakers: micro, cat. no. 13-0051, macro, cat. no. 13-0050
- Type of Sealing : Viton, cat. no. 1000578
- Extraction Thimbles: Type SE33A, 33 x 80 mm, cat. no. 13-0054
Type SE33B, 33 x 94 mm, cat. no. 13-0057
- Extraction Baskets: SHK2, cat. no. 13-0062
- Boiling Stones: cat. no. 1000774
- Air Pressure / Compressor: min. 4.5 bar or cat. no. 13-0010
- Water Pressure or Recirculating Cooler: min. 0.5 bar or FL 601, cat. no. 10-0046 or FL1201, cat. no. 10-0045



CRUDE FAT IN FEED – SPECIAL METHOD FOR FOAMING SAMPLES

| Program Step | Program parameter | Comment |
|------------------------|-------------------|---|
| Safety Temperature | 200 °C | |
| Extraction Temperature | 150 °C | |
| Reduction Interval | 4 min | |
| Reduction Pulse | 2 s | |
| Hot Extraction | 30 min | |
| Evaporation A | 4-5 x Interval | After A the level of solvent should be at least 10 mm below the thimble |
| Rinsing Time | 60 min | |
| Evaporation B | 3-4 x Interval | After B the extraction beaker should be nearly empty |
| Evaporation C | 4 min | |

After the program is run, the extraction beakers are dried in the drying chamber for 1.5 hours at 100°C ± 2 K. Then, they are put into the desiccator, left to cool down to room temperature and weighed with a precision of +/- 1 mg. In order to check the weight consistency, the samples are dried for another 30 minutes and weighed again after cooling down. This procedure is repeated as long until two successive weighings show no more difference than 0.1 % of the initial sample weight. Should the weight increase then the previous lower value should be taken. Extraction, drying, and weighing have to be effected consecutively.

8. Evaluation

8.1. Calculation

w: Crude fat content in g/100 g (equivalent %) of the sample to be calculated as follows:

$$w = \frac{(m_2 - m_1) * 100}{m_0}$$

m1: Weight (g) of the empty extraction beaker with boiling stones

m2: Weight (g) of the extraction beaker with fat after the drying process

m0: Sample weight (g)

The result is given rounded off to one fractional digit.

8.2. Reproducibility

The difference in results of double-determinations should not exceed:

At a crude fat content of:

- less than 5 % - 0.2 % absolute
- 5 % - 10 % - 4.0 % of the higher result
- more than 10 % - 0.4 % absolute

8.3. Comment

Possibly, no hydrolysis is necessary for the crude fat determination in single feedstuffs - this depends on the fat type.



CRUDE FAT IN FEED – SPECIAL METHOD FOR FOAMING SAMPLES

9. Product Table

9.1. Single Feed

| Sample Type | Sample weight [g] +/- 10 % | Crude fat content [%] |
|--|-------------------------------|--------------------------|
| Babssu extraction coarse meal | 10.000 | max. 4 |
| Babssu extraction coarse meal, partly husked | 10.000 | max. 4 |
| Babassu cake | 10.000 | min. 5 |
| Coarse cotton seed extraction meal from husked seed | 10.000 | max. 4 |
| Coarse cotton seed extraction meal from husked seed, fattened | 10.000 | min. 4 |
| Coarse cotton seed extraction meal from partly husked seed treated with formaldehyde, for cattle, sheeps and goats | 10.000 | max. 4 |
| Coarse cotton seed extraction meal from unpeeled seed | 10.000 | max. 4 |
| Coarse cotton seed extraction meal from unpeeled seed, fattened | 10.000 | min. 4 |
| Coarse cotton seed cake from husked seed | 10.000 | min. 5 |
| Coarse cotton seed cake from partly husked seed | 10.000 | min. 5 |
| Ecuador palm kernel cake | 10.000 | max. 4 |
| Peanut, peeled, husked | 10.000 | max. 4 |
| Peanut extraction coarse meal from husked seed | 10.000 | max. 4 |
| Peanut extraction coarse meal from partly husked seed | 10.000 | max. 4 |
| Peanut extraction coarse meal from partly husked seed, fattened | 10.000 | min. 4 |
| Peanut extraction coarse meal from non-husked seed | 10.000 | max. 4 |
| Peanut bran | 10.000 | min. 10 |
| Peanut cake from husked seed | 10.000 | min. 5 |
| Peanut cake from partly husked seed | 10.000 | min. 5 |
| Fish liver meal | 10.000 | max. 10 |
| Fish liver meal, partly hydrolysed | 10.000 | max. 2,5 |
| Fish meal type 55 | 10.000 | max. 12 |
| Fish meal type 60 | 10.000 | max. 12 |
| Fish meal type 64 | 10.000 | max. 12 |
| Meat-and-bone meal | 10.000 | max. 14 |
| Bone coarse | 10.000 | max. 5 |
| Poultry waste, dried | 10.000 | max. 12 |
| Poultry waste, rich in fat, dried | 10.000 | min. 12 |
| Hemp extraction coarse meal | 10.000 | max. 4 |
| Hemp cake | 10.000 | min. 5 |
| Cacao extraction coarse meal | 10.000 | max. 4 |
| Kapok extraction coarse meal | 10.000 | max. 4 |
| Kapok extraction coarse meal, fattened | 10.000 | min. 4 |
| Kapok cake | 10.000 | min. 5 |
| Coconut extraction coarse meal | 10.000 | max. 4 |
| Coconut extraction coarse meal, fattened | 10.000 | min. 4 |
| Coconut cake | 10.000 | min. 5 |
| Copra, dried | 3.500 | min. 60 |
| Lin extraction coarse meal | 10.000 | max. 4 |
| Lin extraction coarse meal, fattened | 10.000 | min. 4 |
| Lin cake | 10.000 | min. 5 |
| Macoya palm kernel extraction coarse meal | 10.000 | max. 4 |
| Macoya palm kernel extraction coarse meal, fattened | 10.000 | min.4 |
| Macoya palm kernel cake | 10.000 | min. 5 |
| Maize germ | 10.000 | min. 20 |
| Coarse maize germ meal (maize milling industry) | 10.000 | max. 4 |
| Coarse maize germ meal (starch industry) | 10.000 | max. 4 |
| Maize germ bran, extracted | 10.000 | max. 4 |

**CRUDE FAT IN FEED – SPECIAL METHOD FOR FOAMING SAMPLES**

| Sample Type | Sample weight [g] +/- 10 % | Crude fat content [%] |
|--|-------------------------------|--------------------------|
| Maize germ cake (maize milling industry) | 10.000 | min. 4 |
| Maize germ cake (starch industry) | 10.000 | min. 4 |
| Almond cake | 10.000 | min. 5 |
| Milk powder, partly skimmed | 10.000 | min.10 |
| Poppy seed extraction coarse meal | 10.000 | max. 5 |
| Poppy seed cake | 10.000 | min. 5 |
| Murumuru Palm kernel extraction coarse meal | 10.000 | max. 4 |
| Niger seed extraction coarse meal | 10.000 | max. 4 |
| Nut cake | 10.000 | min. 4 |
| Olive extraction coarse meal | 10.000 | max. 4 |
| Palm kernel extraction coarse meal | 10.000 | max. 4 |
| Palm kernel extraction coarse meal, fattened | 10.000 | min. 4 |
| Palm kernel cake | 10.000 | min. 5 |
| Rape extraction coarse meal | 10.000 | max. 4 |
| Rape extraction coarse meal, fattened | 10.000 | min. 4 |
| Rape extraction coarse meal with reduced goitrin content | 10.000 | max. 4 |
| Rape extraction coarse meal with reduced goitrin, fattened | 10.000 | min. 4 |
| Rape seed cake | 10.000 | min. 5 |
| Rice germ | 10.000 | min. 17 |
| Rice germ extraction coarse meal | 10.000 | max. 4 |
| Rice germ extraction coarse meal, fattened | 10.000 | min. 4 |
| Rice germ cake | 10.000 | min. 5 |
| Safflower extraction coarse meal from husked seed | 10.000 | max.4 |
| Safflower extraction coarse meal from partly husked seed | 10.000 | max.4 |
| Safflower extraction coarse meal from husked seed | 10.000 | min. 5 |
| Silk worm cocoon coarse meal, extracted | 10.000 | max. 4 |
| Sesame extraction coarse meal | 10.000 | max. 4 |
| Sesame cake | 10.000 | min. 5 |
| Soy extraction coarse meal | 10.000 | max. 4 |
| Soy extraction coarse meal, steamed | 10.000 | max. 4 |
| Soy extraction coarse meal, steamed, fattened | 10.000 | min. 4 |
| Soy cake | 10.000 | min. 5 |
| Sunflower extraction coarse meal from husked seed | 10.000 | max. 4 |
| Sunflower extraction coarse meal from husked seed, fattened | 10.000 | min. 4 |
| Sunflower extraction coarse meal from partly husked seed | 10.000 | max. 4 |
| Sunflower extraction coarse meal from partly husked seed, fattened | 10.000 | min. 4 |
| Sunflower extraction coarse meal from unpeeled seed | 10.000 | max. 4 |
| Sunflower extraction coarse meal from unpeeled seed, fattened | 10.000 | min. 4 |
| Sunflower seed cake from husked seed | 10.000 | min. 5 |
| Sunflower seed cake from partly husked seed | 10.000 | min. 5 |
| Tengkawang coarse meal, extracted | 10.000 | max. 4 |
| Meat-and-bone meal | 10.000 | max. 11 |
| Meat-and-bone meal, rich in fat | 10.000 | min. 11 |
| Tucum extraction coarse meal | 10.000 | max. 4 |
| Uricuri extraction coarse meal | 10.000 | max. 4 |

**CRUDE FAT IN FEED – SPECIAL METHOD FOR FOAMING SAMPLES****7.2. Mixed Feed**

| Sample Type | Sample weight [g] +/- 10 % | Crude fat content [%] |
|--|-------------------------------|--------------------------|
| Milk replacer for growing calves (complete feed)) | 8.000 | 5 – 30 |
| Milk replacer I for fattening calves (complete feed) | 8.000 | 8 – 30 |
| Milk replacer II for fattening calves from approx. 80 kg on (complete feed) | 8.000 | 15 – 30 |
| Supplementary feed rich in energy for skimmed milk for fattening calves | 4.000 | 30 – 60 |
| Dairy concentrate I (supplementary feed for dairy cows) | 10.000 | max. 6 |
| Dairy concentrate II (supplementary feed for dairy cows) | 10.000 | max. 6 |
| Dairy concentrate III (supplementary feed for dairy cows) | 10.000 | max. 8 |
| Dairy concentrate IV (supplementary feed for dairy cows, rich in protein) | 10.000 | max. 8 |
| Cattle feed I (supplementary feed for fattening cattle) | 10.000 | max.6 |
| Cattle feed II (supplementary feed for fattening cattle, rich in protein) | 10.000 | max. 6 |
| Milk replacer for piglet (complete feed) | 10.000 | min 4 |
| Piglet starter I (complete feed) up to approx. 20 kg | 10.000 | max. 7 |
| Piglet starter II (complete feed) up to approx.35 kg | 10.000 | max. 7 |
| Complete feed I for fattening pigs up to approx. 50 kg | 10.000 | max. 8 |
| Complete feed II for fattening pigs from approx. 50 kg on | 10.000 | max. 10 |
| Complete feed for fattening pigs from approx. 35 kg on | 10.000 | max. 9 |
| Complete feed for lactating sows | 10.000 | max. 8 |
| Supplementary feed for suckling piglet | 10.000 | max. 6 |
| Supplementary feed I for fattening pigs | 10.000 | max. 12 |
| Supplementary feed II for fattening pigs | 10.000 | max. 12 |
| Supplementary feed I for stock pigs | 10.000 | max. 12 |
| Milk replacer for ewe lambs (complete feed) | 8.000 | min. 20 |