

Solution Expiration Dates

A GUIDE FOR ENVIRONMENTAL PROFESSIONALS

NOTE: Expiration dates are clearly labeled on the bottle. Expiration times for unopened / unmixed solutions are from date of manufacture, not date of purchase. TruLine solutions are for the TruLine laboratory ISEs.



CONDUCTIVITY

Unopened - Glass Quarts	1 Year
Unopened - Plastic Pints	18 Months
Opened	1 Month

TruLine ISE CALIBRATION, ISA & FILLING SOLUTIONS

Unopened	15 Months
Opened	6 Months

TURBIDITY

Unopened	1 Year
Opened	6 Months

ZOBELL

Unmixed	5 Years
Mixed	6 Months

AMMONIUM

Unopened	18 Months
Opened	6 Months

pH

Unopened	2 Years
Opened	2 Months

NITRATE

Unopened	1 Year
Opened	6 Months

CONFIDENCE SOLUTION

Unopened	2 Years
Opened	3 Months

O₂ PROBE ELECTROLYTE SOLUTION

Unmixed Crystals	No Stated Date
Mixed Solution	1 Year

YSI, a Xylem brand
1725 Brannum Lane
Yellow Springs, OH 45387

+1.937.767.7241
info@ysi.com
YSI.com

Who's
Minding
the Planet?®



a xylem brand