



Who's
Minding
the Planet?®

SOLUTION EXPIRATION DATES

CONDUCTIVITY

Unopened -- Glass Quarts _____ 1 Year
Unopened -- Plastic Pints _____ 18 Months
Opened _____ 1 Month

ZOBELL

Unmixed _____ 5 Years
Mixed _____ 6 Months

pH

Unopened _____ 2 Years
Opened _____ 6 Months

CONFIDENCE SOLUTION

Unopened _____ 2 Years
Opened _____ 3 Months

TURBIDITY

Unopened _____ 1 Year
Opened _____ 6 Months

AMMONIUM

Unopened _____ 18 Months
Opened _____ 6 Months

NITRATE

Unopened _____ 1 Year
Opened _____ 6 Months

O₂ PROBE ELECTROLYTE SOLUTION

Unmixed Crystals _____ No Stated Expiration
Mixed Solution _____ 1 Year

TruLine ISE CALIBRATION, ISA, AND FILLING SOLUTIONS

Unopened _____ 15 Months
Opened _____ 6 Months

NOTE:

Expiration times for unopened / unmixed solutions are from date of manufacture, not date of purchase.

Expiration dates are clearly labeled on the bottle.

TruLine solutions are for the TruLine laboratory ISEs.

