



ProSolo Ordering Guide

WATER QUALITY FIELD INSTRUMENT

Step 1: Select Handheld

All handhelds include a rechargeable lithium-ion battery, hand strap, USB cable, universal AC charger, data export cable, Quick Start Guide, and USB memory stick containing KorDSS software and digital copy of the User Manual.

- **626650**: ProSolo Handheld



Key Features:

- Color display
- 100,000+ data set memory
- Data export to USB flash drive
- Rechargeable battery
- Kor Software

ProSolo Handheld

Step 2: Select Cable Assembly

Integrated Probe and Cable Assemblies

Simplify instrument setup with integrated probe and cable assemblies. Available in a variety of parameter-specific options.

- **627150**: ODO/CT cable with integral optical dissolved oxygen (with pre-installed sensor cap), conductivity, and temperature sensors; **1, 4, 10, 20, 30, 50,** or **100** meter length
- **627200**: ODO/T cable with integral optical dissolved oxygen (with pre-installed sensor cap) and temperature sensors; **1, 4, 10, 20, 30, 50,** or **100** meter length
- **626400**: ProOBOD self-stirring optical BOD probe/cable assembly, lab probe; 115 V; **3** meter length
- **626401**: ProOBOD self-stirring optical BOD probe/cable assembly, lab probe; 230 V; **3** meter length



ODO/CT



ODO/T



OBOD

Step 3: Select Standards

- **065270**: Conductivity standard, 1000 umhos/cm (quart)
- **065272**: Conductivity standard, 10000 umhos/cm (quart)
- **065274**: Conductivity standard, 100000 umhos/cm (quart)
- **060907**: Conductivity calibrator, 1000 umhos/cm (8 ea, pint)
- **060906**: Conductivity calibrator, 1413 umhos/cm (8 ea, pint)
- **060911**: Conductivity calibrator, 10000 umhos/cm (8 ea, pint)
- **060660**: Conductivity calibrator, 50000 umhos/cm (8 ea, pint)



065270



065272



065274

Step 4: Select Accessories

- 626946:** Large, hard-sided carrying case (Fits 10, 20, and 30 meter cables)
- 603075:** Large, soft-sided carrying case
- 626945:** Small, hard-sided carrying case (Fits 1 and 4 meter cables plus a flow cell)
- 603076:** Flow cell for ODO/CT cables (requires single port adapter; 603078)
- 603078:** Adapter required for ODO/CT flow cell (603076)
- 603056:** Flow cell mounting spike
- 063507:** Tripod
- 063517:** Ultra clamp
- 603070:** Shoulder strap
- 603069:** Belt clip
- 626918:** 1 lb weight - stackable on guard
- 605978:** 4.9 oz weight - stackable on guard
- 626942:** USB car charger
- 626943:** Small external USB portable battery pack
- 626944:** Large external USB portable battery pack



Smart, Hard-Sided Carrying Case



Stackable Guard Weights



Tripod

Replacement Parts

- 626846:** Rechargeable Lithium-ion battery pack
- 626969:** USB flash drive with user manuals and KorDSS
- 626991:** USB cable for charging and PC connection
- 626992:** Cable for connection to USB stick
- 627195:** Calibration cup for ODO/CT cables
- 603062:** Cable management kit
- 626940:** AC charger (USA). Includes power supply and USB cable (626991)
- 626941:** AC charger (International). Includes power supply, USB cable and outlet adapters
- 627180:** ODO Sensor cap for ODO/CT and ODO/T assemblies
- 626482:** ODO Sensor cap for ProOBOD assemblies (626400 or 626401)



Interested in ordering a system but still have questions?

Please use the contact information below. We will ensure the system you order is right for your application.



Simplify your setup with the ProSolo and an integrated probe and cable assembly.



YSI, a Xylem brand
1725 Brannum Lane
Yellow Springs, OH 45387

- +1.937.767.7241
- info@ysi.com
- YSI.com

© 2020 Xylem, Inc. XA00066 0220



[YSI.com/ProSolo](https://www.ysi.com/ProSolo)