


# Kor Software

VERSION 1.3.5.0

**Kor Software** is an accessible, feature-rich platform for managing YSI instrumentation and water quality data. Easily calibrate sensors, setup long-term deployments, and ensure equipment is operating at peak performance.



### Compatible Instruments:

Kor Software replaces both KorEXO and KorDSS as the primary user interface for managing EXO and Pro Series instruments and data.

- EXO Sondes
- EXO Handheld
- ProSwap Logger
- ProDSS Handheld
- ProSwap Handheld
- ProSolo Handheld

### Kor Software – v1.3.5.0

November 2023

### What's New:

Bug fixes and added export functionality.

### Change Notes:

- Bug fixes for database imports
- Bug fixes for handheld downloads
- Improvements to Modbus enable/disable interface for ProSwap Logger
- Added export function for Pro Series instrument configuration
- Added print and export function for EXO NitraLED site correction and calibration
- User Interface corrections

## PREVIOUS UPDATES

### Kor Software – v1.1.8.0

September 2021

#### What's New:

This release corrects an issue with the ProDSS database import that previously resulted in a crash.

#### Change Notes:

- Fixed issue in DSS database import process
- Added routine to existing databases

### Kor Software – v1.1.6.0

July 2021

#### What's New:

Kor Software replaces both KorEXO and KorDSS as the primary user interface for managing EXO and ProDIGITAL (ProSwap, ProSolo, and ProDSS) instruments and data.

#### Change Notes:

- Allow downloading of partial deployment file
- Now displaying ISE millivolts during calibration process and in calibration reports
- Fixed crash caused by duplicated site names
- Fixed crash during ProDSS Handheld data download
- Fixed crash due to import of old wiper calibration
- Fixed crash related to ProSolo Handheld connection
- Fixed crash when creating configuration template from ProDSS Handheld
- Removed "Missing CT" message when calibrating ODO/CT sensor
- Fixed issue where download could not be cancelled
- Various typo and other bug fixes

YSI, a Xylem brand  
1725 Brannum Lane  
Yellow Springs, OH 45387

+1.937.767.7241  
info@ysi.com  
YSI.com

Who's  
Minding  
the Planet?®



a xylem brand